



旋转痰盂设计，方便患者使用，豪华型陶瓷痰盂缸可拆卸，便于消毒清洗
Luxurious rotating ceramic design, convenient for the patient to use and easy to take down to sterilize and clean



2个水瓶，其中一个水瓶装消毒水，手机管线消毒，防止病人交叉感染。
2 water bottles, one of which is equipped with sterilized water, sterilized by dental handpiece tubes to prevent cross infection.



多功能脚踏开关
Multi function foot control, can directly control the spittoon



反射式 LED 口腔灯
采用先进的 LED 芯片技术比传统卤素灯照度更高更接近自然光，可感应或手动开关控制强弱光
Reflective LED oral light, using advanced LED chip technology, higher than traditional halogen lighting, closer to natural light, induction or manual switch to control light and weak light

DK-18 牙科综合治疗台 INTEGRAL DENTAL UNIT

Specifications:

Power Voltage: 220V, 50HZ
Motor Voltage: 24V
Water Pressure: 0.2-0.4Mpa
Air Pressure: 0.5-0.8Mpa

电流电压：220V 50HZ
电机电压：24V
水源水压：0.2MPa-0.4MPa
气源气压：0.5MPa-0.8MPa

Rotable glass/ceramic spittoon	1set
Auto spittoon flush and cup filler control system	1set
24V noiseless dc motor chair	1set
3-Way syringe (cold /hot water)	2set
Led operating light	1set
Saliva ejector with high and low suction	1set
Dental handpiece tube disinfection system	1set
Water tank (air pressurized) (built-in)	2set
Instrument tray with air brake	1set
Operating stool	1set

可旋转可拆卸陶瓷或玻璃痰盂	1套
冲盂漱口定量给水自动控制系统	1套
24V 静音直流电动椅	1套
3用喷枪（冷/热）	2支
LED 口腔冷光灯	1套
强弱吸唾系统	1套
内置式净水系统 消毒水系统	2套
气刹装置工具盘	1套
医生转椅	1张





MYD-108 牙科综合治疗台

INTEGRAL DENTAL UNIT

Specifications:

Power Voltage:220V,50HZ
Motor Voltage:24V
Water Pressure:0.2-0.4Mpa
Air Pressure:0.5-0.8Mpa

Rotable ceramicspittoon	1set
Auto spittoon flush and cup filler control system	1set
24V noiseless dc motor chair	1set
3-Way syringe	1set
Operating light (optical led light)	1set
Saliva ejector with high and low suction	1set
Water tank(air pressurized)(built-in)	1set
Instrument tray with air brake	1set
Operating stool	1set

电流电压：220V 50HZ
电机电压：24V
水源水压：0.2MPA-0.4MPA
气源气压：0.5MPA-0.8MPA

可旋转陶瓷痰盂	1套
冲盂漱口定量给水自动控制系统	1套
24V 静音直流电动椅	1套
3用喷枪	1支
口腔冷光灯	1套
强弱吸唾系统	1套
内置式净水系统 消毒水系统	1套
气刹装置工具盘	1套
医生转椅	1张



正面展示图
FRONT VIEW



背面展示图
BACK VIEW



可任意搭配多款 LED 灯
A variety of LED lights for your choice



配备透明保护罩
Protective cover



标配副控展示，可加装 3 用枪以及副控按键
Standard control, it can be installed with three way syringe and electronic push-button



旋转痰盂设计，方便患者使用
Luxurious rotating ceramic design, convenient for the patient to use.



多功能脚踏开关
Multi function foot control .Better durable



MYD-1000

MYD-1000 标准型牙科椅介绍

DENTAL CHAIR UNIT


可旋转可拆卸陶瓷或玻璃痰盂	1 套
冲盂漱口定量给水自动控制系统	1 套
24V 静音直流电动椅	1 套
3 用喷枪 (冷 / 热)	2 支
LED 口腔冷光灯	1 套
强弱吸唾系统	1 套
内置式净水系统	1 套
气刹装置工具盘	1 套
医生转椅	1 张

90 degree rotatable frosted glass spittoon	1set
Auto spittoon flush and cup filler control system	1set
24V noiseless DC motor chair	1set
3-way syringe (cold /hot water)	2set
LED operating light	1set
Saliva ejector with high and low suction	1set
Water tank(air pressurized)(built-out)	1set
Instrument tray with air brake	1set

Specifications:

Power Voltage:220V,50HZ
Motor Voltage:24V
Water Pressure:0.2-0.4Mpa
Air Pressure:0.5-0.8Mpa

电流电压 : 220V 50HZ
电机电压 : 24V
水源水压 : 0.2MPa-0.4MPa
气源气压 : 0.5MPa-0.8MPa



MYD-208

Configuration Introduced

Ledfilm viewer	1set
24V noiseless dc motor	1set
Instrument tray with air brake	1set
Led colo light lamp	1set
Water tank(air pressurized)(build-outside)	1set
Automatically flushing spittoon control system	1set
Built-function foot control	1set
Saliva ejector with high and low suction	1set
3 Way syringe(hot/cold each)	1set each

For option

Im ported or china original scaler
Im ported or china original curing light
Im ported or china original 4 holes high-speed handpiece
Im ported or china original 4 holes low-speed handpiece
Led bleaching system
Air compressor

Technology parameters

Input power:220v/50hz
Motor voltage:24v
Water pressure:0.2Mpa-0.4Mpa
Air pressure:0.5Mpa-0.8Mpa



MYD-308

Configuration Introduced

Luxury electrical dental chair(with compensation function)	1set
24V noiseless dc motor	1set
Rotable tempered glass spittoon	1set
Led colo light lamp	1set
Computer control system with 3 memory set	1set
Automatically flushing spittoon control system	1set
Multi-function foot control	1set
Strong/weak suction with filter	1set
3 Way syringe(hot/cold each)	1set each

For option

Im ported or china original scaler
Im ported or china original curing light
Im ported or china original 4 holes high-speed handpiece
Im ported or china original 4 holes low-speed handpiece
Led bleaching system
Air compressor

Technology parameters

Input power:220v/50hz
Motor voltage:24v
Water pressure:0.2Mpa-0.4Mpa
Air pressure:0.5Mpa-0.8Mpa



MYD-008

Configuration Introduced

3-way syringe(cold &hot)	1set each
high suction and saliva ejector	1set each
operation light	1set
ceramic cuspidor	1set
LED film viewer	1set
automatic bowl rinsing and cup filling	1set
water purified system of hand pieces	1set
assistant stand	1set
foot switch	1set
doctor's chair	1set
PU cushion	1set

For option

Can change for top mounted
Oilless air compressor
High speed handpiece
Low speed handpiece
Digital intra-oral camera
LED dental curing light
Polisher
Scaler

Technology parameters

Common Technical Parameters Introduced
Power Supply AC230V 50Hz(customizable)
Motor Power DC24V
Air Pressure 0.5MPa - 0.8MPa
Water Pressure 0.2MPa - 0.4MPa



MYD-C01



MYD-C02

一拖一空压机 ONE TO ONE AIR COMPRESSOR

技术参数
 额定电压：220VAC 50HZ
 电流：2.4A
 功率：545W
 容积流量：70L/min
 启动压力：0.5Mpa
 额定排气压力：0.8Mpa
 储气罐容积：32L
 噪音：55-62dB
 重量：28kg
 包装尺寸：46x46x69cm(纸箱包装)

Features
 Voltage: 220VAC 50Hz
 Electric current: 2.4A
 Power: 545W
 Frequency: 70L/min
 Start pressure: 0.5Mpa
 Capacity pressure: 0.8Mpa
 Air tank volume: 32L
 Noise: 55-62dB
 G.W.: 28kg
 Packing size: 46x46x69cm (carton)



MYD-C05



MYD-C06

一拖三空压机 ONE TO THREE AIR COMPRESSOR

技术参数
 额定电压：220VAC 50HZ
 电流：4.8A
 功率：1090W
 容积流量：140L/min
 启动压力：0.5Mpa
 额定排气压力：0.8Mpa
 储气罐容积：60L
 噪音：65-68dB
 重量：66kg
 包装尺寸：90x44x77cm(木箱包装)

Features
 Voltage: 220VAC 50Hz
 Electric current: 4.8A
 Power: 1090W
 Frequency: 140L/min
 Start pressure: 0.5Mpa
 Capacity pressure: 0.8Mpa
 Air tank volume: 60L
 Noise: 65-68dB
 G.W.: 66kg
 Packing size: 90x44x77cm (wooden box)

一拖二空压机 ONE TO TWO AIR COMPRESSOR

技术参数
 额定电压：220VAC 50HZ
 电流：3.7A
 功率：840W
 容积流量：105L/min
 启动压力：0.5Mpa
 额定排气压力：0.8Mpa
 储气罐容积：32L
 噪音：58-64dB
 重量：30kg
 包装尺寸：46x46x71cm(纸箱包装)

Features
 Voltage: 220VAC 50Hz
 Electric current: 3.7A
 Power: 840W
 Frequency: 105L/min
 Start pressure: 0.5Mpa
 Capacity pressure: 0.8Mpa
 Air tank volume: 32L
 Noise: 58-64dB
 G.W.: 29kg
 Packing size: 43x43x70cm (carton)



MYD-C03



MYD-C04

一拖四空压机 ONE TO FOUR AIR COMPRESSOR

技术参数
 额定电压：220VAC 50HZ
 电流：7.4A
 功率：1680W
 容积流量：210L/min
 启动压力：0.5Mpa
 额定排气压力：0.8Mpa
 储气罐容积：60L
 噪音：65-68dB
 重量：70kg
 包装尺寸：90x44x77cm(木箱包装)

Features
 Voltage: 220VAC 50Hz
 Electric current: 7.4A
 Power: 1680W
 Frequency: 210L/min
 Start pressure: 0.5Mpa
 Capacity pressure: 0.8Mpa
 Air tank volume: 60L
 Noise: 65-68dB
 G.W.: 70kg
 Packing size: 90x44x77cm (wooden box)



MYD-C07



MYD-C08



MYD-C09

一拖五空压机

ONE TO FIVE AIR COMPRESSOR

技术参数

额定电压：220VAC 50HZ
电流：7.2A
功率：1635W
容积流量：210L/min
启动压力：0.5Mpa
额定排气压力：0.8Mpa
储气罐容积：90L/128L
噪音：68-72dB
重量：89kg
包装尺寸：105x47x78cm(木箱包装)

Features

Voltage: 220VAC 50Hz
Electric current: 7.2A
Power: 1635W
Frequency: 210L/min
Start pressure: 0.5Mpa
Capacity pressure: 0.8Mpa
Air tank volume: 90L
Noise: 68-72dB
G.W.: 89kg
Packing size: 105x47x78cm (wooden box)



MYD-C10

一拖六空压机

ONE TO SIX AIR COMPRESSOR

技术参数

额定电压：220VAC 50HZ
电流：11A
功率：2520W
容积流量：315L/min
启动压力：0.5Mpa
额定排气压力：0.8Mpa
储气罐容积：90L/128L
噪音：72-75dB
重量：95kg
包装尺寸：105x47x78cm(木箱包装)

Features

Voltage: 220VAC 50Hz
Electric current: 11A
Power: 2520W
Frequency: 315L/min
Start pressure: 0.5Mpa
Capacity pressure: 0.8Mpa
Air tank volume: 90L
Noise: 72-75dB
G.W.: 95kg
Packing size: 105x47x78cm (wooden box)



MYD-C11

负压泵 (一拖二)

SUCTION UNIT (ONE TO TWO)

可带治疗机数：2 ~ 3 台
技术参数 SVS300:
电源电压 (VAC)：220V
频率 (Hz)：50/60
功率 (Kw)：0.7/0.8
电流 (A)：4.8/4.1
真空度 (m bar)：150/160
电机转速 (rpm)：2,920/3,450
可带治疗机数：2-3 (同时使用 3 台治疗机可达 60%的效率)
重量 (kg)：20/25 (连包装)

Usage: 2 ~ 3 dental units
Specification:
Voltage: 220V AC 50/60Hz
Power(Kw): 0.7/0.8
Electric current: 4.1/4.8
Vacuum rate(m bar): 150/160
Rotation speed(rpm): 2,920/3,450
G.W.:25kg
Dimensions(LxWxH):45x50x67cm

选位阀 1

选位阀 2



MYD-C12

便携负压系统

DENTAL SUCTION

1. With an air/liquid separation unit, it used to separates the liquid and solid particles
2. with drain pump system, it discharges the residue of liquid and blood automatically after stopping the vacuum pump/motor. it prevent clogging problems

Advantages

1. With display crystal screen to show the liquid position of the container, and other working state.
2. Choose drain water by automatic or manual or manual from buzzer/ channel
3. This suction offer one dental chair.

Features

1. Powerful Source
2. The Second Protector For Dentist
3. Small Space, High Tech
4. Prevention Of Infection
5. Pressure under the control of the extraordinary power.
6. Reliable Technology
7. Simple Maintain
8. Flexibility

Technical Parameter

power 110VAC/60HZ;
220VAC/50HZOption
Inputpower 370W
Maximumwinpressure 11.5kPa
Vacuumdegree -11kPa
Weight 26.5kg
Noise 53dB
Motorrevolution 2800r/min



MYD-C13
拖一泵头 (545w)
One to one motor



MYD-C14
拖二泵头 (840w)
One to two motor



MYD-C15
泵头外壳
Motor cover



MYD-C16
泵头外壳 (带出线口)
Motor cover



MYD-C17
左缸盖
Cylinder cover



MYD-C18
偏心轮 (拖一)
Wheel



MYD-C19
偏心轮 (拖二)
Wheel



MYD-C20
转子轴承
Rotor bearing



MYD-C21
偏心轮轴承 (拖一)
Wheel bearing



MYD-C22
偏心轮轴承 (拖二)
Wheelbearing



MYD-C23
连杆 (拖一)
Connecting rod



MYD-C24
连杆 (拖二)
Connecting rod



MYD-C25
皮碗压板 (拖一)
Cup plate



MYD-C26
皮碗压板 (拖二)
Cup plate



MYD-C27
缸套 (拖一)
Cylinder liner



MYD-C28
缸套 (拖二)
Cylinder liner



MYD-C29
阀板
Valve plate



MYD-C30
皮碗
Seal cup



MYD-C31
气门阀片
Air lock



MYD-C32
阀片固定板
Lock plate



MYD-C33
阀片固定板
Lock plate



MYD-C34
泵头防尘罩
Dust Protect cover



MYD-C35
泵头冷却风扇
Cooling fans



MYD-C36
阀板密封圈
Valve Plate seal ring



MYD-C37
缸盖密封圈 I
Cylinder



MYD-C38
消音器
Noise reduction



MYD-C39
空气滤芯
Air filter



MYD-C40
消音器
Noise reduction



MYD-C41
立式压力开关
Pressure switch



MYD-C42
卧式压力开关 (a)
Pressure switch



MYD-C43
卧式压力开关 (b)
Pressure switch



MYD-C44
压力表
Pressure gauge



MYD-C45
排气电磁阀
Solenoid valve



MYD-C46
排气电磁阀
Solenoid valve



MYD-C47
球阀开关 (a)
Switch



MYD-C48
球阀开关 (b)
Switch



MYD-C49
球阀开关 (c)
Switch



MYD-C50
单向阀 (b)
One way valve



MYD-C51
单向阀 (c)
One way valve



MYD-C52
单向阀 (d)
One way valve



MYD-C53
单向阀 (e)
One way valve



MYD-C54
安全阀
Safety valve



MYD-C55
排水阀
Water valve



MYD-C56
高压进气直管
Intake pipe



MYD-001

MYD-002

便携椅介绍

PORTABLE DENTAL CHAIR (LED LIGHT AND FOLDABLE)

产品具有轻便灵活，可任意方向移动定位，结构牢固稳定，安装拆卸简便，采用优质结构钢喷漆工艺制造。产品座椅面可调高低位，从离地400-500mm，靠背可调后倾角度从105-160度。承载为135kg。

The portable chair has the advantages of portability and flexibility, and can be moved and positioned in any direction. The structure is firm and stable, and the installation and disassembly are simple and convenient, and the structure is made of high quality structural steel spraying technology product seat surface can be adjusted high and low, from the ground 400-500mm, backrest adjustable, tilt angle from 105-160 degrees. Load is 135kg



MYD-003

挂墙式涡轮机 (a)
Metal Turbines (wall type)



MYD-004

涡轮机 (b)
Metal Turbines



MYD-005

涡轮机 (c)
Plastic Turbines



MYD-006

涡轮机 (d)
Plastic Turbines with 2HP



MYD-007

便携式涡轮机 拉杆箱 内置储气罐
PORTABLE DENTAL UNIT BUILT-IN STORAGE TANKS

Standard Accessories

Built-in air compressor	1set
Built-in storage tanks	1set
Saliva ejector	1set
3-way air/water syringe	1set
High speed air turbine Handpiece tube	1set
Low speed air turbine Handpiece tube	1set
Foot pedal	1set
Clean water bottle 600ML	1set

Specifications:

Power: 550W
Net /Gross Weight: 21Kg/23Kg
Dimension/CM: L50*W27*H45
Package size/CM: L53*W35*H53

Options

Built-in LED curing light
Built-in ultrasonic scaler-Woodpecker
High speed handpiece
Low speed handpiece



MYD-008

便携式涡轮机 拉杆箱
PORTABLE DENTAL UNIT

Standard Accessories

Built-in air compressor	1set
Saliva ejector	1set
3-way air/water syringe	1set
High speed air turbine Handpiece tube	1set
Low speed air turbine Handpiece tube	1set
Foot pedal	1set
Clean water bottle 600ML	1set
3-way syringe bottle	1set
Waste water bottle	1set

Specifications:

Power: 550W
Net /Gross Weight: 23.5 Kg/?Kg
Dimension/CM: L 43*W34*H 64

Options

Built-in LED curing light
Built-in ultrasonic scaler-Woodpecker
High speed hanpiece
Air motor set

This device can use 【bottle water】 and 【tap water】。
Product dimension: 58x44x84cm (height can be adjusted to 9cm)

Standard Accessories

2 holes/4 holes hand piece tube	2set
Pneumatic weakly absorbing	1set
Triple-purpose gun	1set
600ML water tank	1set
Round foot switch	1set
8 øwater tubing (1.6m)	2set
12 ødewatering soft tubing (1.2m)	1set

Specifications:

Net weight: 18.5kg
Gross weight: 23kg
Package size: 69x52x41 cm (Honeycomb carton)

Options

High speed hand piece
Low speed hand piece
Tooth cleaning machine
Curing light machine
Strong absorption



MYD-009

移动式治疗机 / 小推车 (不配气泵)
MOBILE DENTAL UNIT / TROLLEY (NO COMPRESSOR)



MYD-010

Automatic sterilizer

European class b. The utility model is suitable for disinfection and sterilization of various instruments. Led display, easy to operate, plastic open water tank, high-density bacterial filter. Automatic sterilization, automatic drying, thermal insulation, accurate temperature control, safety protection design and monitoring system

Capacity	12l	18l	23l
Pot size	Φ20*36cm	Φ25*35cm	Φ25*45cm
Display	LED	LED	LED
Product size	66*44.5*39	66*44.5*39	70*44.5*39
Net /gross weight	53/45kgs	55/47kgs	58/50kgs
Power	220v50hz 1.7Kw	220v50hz 1.7Kw	220v50hz 1.7Kw
Output	Printer (optional)	Printer (optional)	Printer (optional)



MYD-015
17 寸多媒体高清一体机
17-Inch multimedia HD



MYD-016
17 寸多媒体高清一体机
17-Inch multimedia HD



MYD-011

Automatic sterilizer

European CLASS B. No water tank design, direct access to water sources, wastewater discharged directly to the sewer, can directly use RO water or pure water. LED display, easy to operate, steam generated in the pot, automatic sterilization, automatic drying, heat preservation function, sterilization time adjustable, precise temperature control, safety protection design, easy maintenance

Capacity	12l	18l	23l
Pot size	Φ20*36cm	Φ25*35cm	Φ25*45cm
Display	LED	LED	LED
Product size	66*44.5*39	66*44.5*39	70*44.5*39
Net /gross weight	45/37kgs	49/41kgs	52/43kgs
Power	220v50hz 1.7Kw	220v50hz 1.7Kw	220v50hz 1.7Kw
Output	Printer (optional)	Printer (optional)	Printer (optional)

技术参数

图像传感器：索尼 1/4 CCD 传感器
系统方案：嵌入式操作系统
显示画面：单画面 / 四画面
系统操作方式：红外线遥控 / 按键控制
备份方式：自动备份到 U 盘
储存容量：最大支持 32G 存储
图片存储格式：JPG
手拍聚焦范围：5MM-50MM
功耗：30W
可选 WIFI 传输：直接和电脑无线连接
高分辨率高低照度，高感光度高质量自动对焦功能
6 颗高亮光长寿命白光 LED
17 寸白色显示器 高档大气

Specifications:

- 1.The SONY 1/4 image sensor CCD sensor
- 2.The system:embedded operating system .
3. The display screen:single picture / four picture
- 4.The operation system:infrared remote control/key control
5. The backup mode:automatic backup to U disk
- 6.The storage capacity:the maximum support 32G memory
- 7.The image storage format:JPG
- 8.The hand focusing range:5MM-50MM
- 9.The power:30W
- 10.The WIFI connection and the computer direct wireless transmission
- 11.Hi-lite LED,6pcs long working life LED,see more detail of the inner mouth clearly
12. 17inch white monitor rich atmosphere
13. High resolution, high and low illumination, high sensitivity, high quality auto focusing function



MYD-012 蒸馏水机
Water distiller

SPECIFICATIONS:

- 1.Voltage: AC110-120V 50-60Hz, AC220-240V 50-60Hz
- 2.Heater efficiency: >95%
- 3.Water measurement: 1L/H
- 4.Capacity: 4L
- 5.Size: 290*260*450mm
- 6.G. W: 3.5kg



MYD-014
封口机
Sealing machine



MYD-013 封口机
Sealing machine

SPECIFICATIONS:

- 1.Voltage: 220/50 V/Hz
- 2.Power: 500W
- 3.Seal width: 250mm(Max)
- 4.Size: 80x410x290mm
- 5.N.W: 6.5KG
- 6.G.W: 8.5KG

- 5.5CM*200M
- 7.5CM*200M
- 10CM*200M
- 15CM*200M
- 20CM*200M
- 25CM*200M




MYD-017
USB Camera



MYD-018
TV Camera



MYD-019
USB Camera



MYD-020
豪华内窥镜支架
Luxury monitor holder



MYD-021
塑料内窥镜支架
Plastic monitor holder



MYD-022
标准内窥镜支架
Monitor holder(METAL)



MYD-023
简易内窥镜支架
Monitor holder(METAL)

Features:

1. Use the international advanced technology, high-efficiency clear, imaging low leak radial, it is only 1% of national standard.
2. Light touched button, controlled with microcomputer, only press one button to select the expose parameter fast and accurately
3. Developing in light room, image in one minute, the doctor can diagnose conveniently for maximum
4. Pneumatic up-down seat, more convenient and more comfortable
5. It is for the use of the digital imaging system mouth

Technical Specification:

1. Power supply: AC220V \pm 10% 50/60HZ 1KVA
2. Tube voltage: 70KV
3. Tube current: 8mA
4. Focus size: 0.8mm
5. Total filtration: 2.5mmAl
6. Exposure time: 0.1-4sec
7. Leak radiation: outside one meter <0.002mGy/h



MYD-024
移动式 X 光机
Mobile X ray unit

MYD-025
防辐射服
Radiation-proof clothes



MYD-026
便携式 X 光机 (a)
Portable X ray unit



MYD-027
便携式 X 光机 (b)
Portable X ray unit



MYD-028
便携式 X 光机 (c)
Portable X ray unit



MYD-029
电子传感器
Digital sensor



MYD-030
传感器夹
Digital sensor holder



MYD-031
胶片定位器
Film locator



MYD-032
明室 X 光片
Bright room X ray film



MYD-033
明室 X 光片
Bright room X ray film



MYD-034
Light curing unit

MYD-035
Built-in light curing

Features:

1. Cover by green plastic shell
2. 5W big power led
3. High articulation optical fiber
4. Working modes: strong. Gradually strong. Flashing
5. Can set the working time, reverse time
6. Bulky battery
7. Solidify all brands of resin materials
8. Also be plugged in power and worked directly

Technology data:

1. Capability of battery: 2200ma/h 3.7V
2. Working voltage: ac100-240v 50-60hz
3. Wavelength: 450-480nm
4. Light power: more than 1000mw/cm²
5. Solidification time and depth: 10s/3mm



Curing light tips (10mm)

Curing light tips (12mm)

Curing light tips (13.8mm)

MYD-036
LED light curing unit

Technical:

1. Unique elbow design
2. original imported high power LED, two types of power 5W, 7W to choose
3. Output power adjustment: 50%/70%/90%/100%
4. Digital LED display, four kinds of time: 10, 20, 30, 40S to choose
5. Brightness 5W \geq 1400mW/cm² 7W \geq 1700mW/cm²
6. Three types of power to power to choose: A type: Li-battery power, B type: AC-110-240V power, C type: AC24V POWER
7. A variety of alternative work patterns: A strong light B gradually bright C pulse
8. Automatic shut-down, low-power warning
9. Sleep function, saving battery



MYD-037
卤素灯泡
BULB



MYD-038
Halogen light curing unit



MYD-039
LED light curing



MYD-040



MYD-041



MYD-042



MYD-043



MYD-044



MYD-045
根管治疗仪带根测功能
Root canal therapy apparatus with apex locator

1. 超大液晶显示屏，清晰读数
2. 设有 10 种程序（设定记忆）供使用
3. 可以根据不同根管扩大针的扭矩设定不同的扭矩值，实时扭矩以条码显示，一目了然，方便医生操作。
4. 在自动反转功能开启状态下，当根管扩大针超过预先设定的扭矩值时，根管扩大针自动反转并退出根管，阻力去除后再自动正转。
5. 治疗手柄由脚踏开关控制。
6. 治疗手柄的马达设计为软启动和软停止，为此没有冲击感，并保证根管扩大针瞬间冲击时不断针。
7. 为保持电池的使用寿命及使用时间，本产品持续 8 分钟未操作时，则自动关闭电源。
8. 弯机头及根管扩大针支架和口角挂钩可承受 135°C 高温高压灭菌消毒。
9. 内置充电电池

1. Large lcd screen, clear reading
2. There are 10 programs (set memory) for use
3. According to different root canal to expand the torque of the needle, set different torque values, real-time torque to bar code display, at a glance, convenient for doctors to operate.
4. In the automatic reversal function open, when the root canal expansion needle exceeds the preset torque value, the root canal expansion needle automatically reverses and exits the root canal, the resistance is removed and then automatically transferred.
5. The treatment handle is controlled by foot switch.
6. The motor of the treatment handle is designed for soft start and soft stop, so there is no shock feeling, and the root canal is expanded. The needle is constantly shot when it is momentarily impacted.
7. In order to maintain the battery life and use time, this product lasts 8 minutes, when not operating, then automatically shut down the power supply. The utility model relates to a 8. Bending head and a root canal expanding needle holder and a mouth hook, which can withstand sterilization at 135 deg c under high temperature and high pressure.
9. Built in rechargeable battery



MYD-046
无线型根管治疗仪
Wireless root canal therapy apparatus

1. 无线手持设计，简单方便
2. 液晶屏幕操作简单
3. 原装瑞士电动马达动力强劲
4. 大容量可充电电池
5. 1:1 弯机
6. 多种反转模式
7. 自动关机功能：8 分钟未使用自动关机
8. 转速扭矩可调
9. 弯机头及根管扩大针支架和口角挂钩可承受 135°C 高温高压灭菌消毒

技术参数
电池电压：dc3.7V 11. 电池容量：1500mah
电压范围：ac110-240v 50/60hz
输出：dc5v 1 a

1. Wireless handheld design, easy and convenient
2. Lcd screen, simple to use it.
3. Original swiss electric motor, power strong
4. High-capacity rechargeable battery
5. 1:1 Contral angel
6. A variety of antiport
7. Automatic shutdown: 8 minutes without use, automatic shutdown. Power saving
8. Adjusted speed, torque
9. The utility model relates to a bent head and a root canal expanding needle holder and a mouth hook, which can withstand sterilization at 135 deg c under high temperature and high pressure

Technical parameter
Battery voltage：dc3.7V
Battery capacity：1500mah
Voltage range：ac110-240v 50/60hz output:dc5v 1a



MYD-047
根管长度测量仪
Apex locator

功能特点：
强大的抗干扰功能：不论根管内有无残髓、血、脓、生理盐水、H2O2 都能测量精确。

显示屏：4"LCD 液晶屏
电源：3 节 AA 电池
功耗：<0.4W
尺寸：95mm*120mm*36mm
重量：175g

Feature:
Strong anti-interference function: no matter the root canal have residual pulp, blood pus, physiological saline H2O2, it can be measured accurately

Screen: 4" LCD
Power: 3 Pieces AA batteries
Power consumption: <0.4W
Size: 95mm*120mm*36mm
Weight: 175g



MYD-048
Color apex locator



MYD-049
White apex locator



MYD-050
Color apex locator



MYD-051
牙龈尖切割器
Gutta cutter



MYD-052
牙髓活力测量仪
Pulp tester



MYD-053
夹持器
Holder



MYD-054
唇钩
Lip hook



MYD-055
注油机
Cleaning lubricator



MYD-056 1000ML Handpiece Oil
MYD-057 500ML Handpiece Oil
MYD-058 250ML Handpiece Oil
MYD-059 100ML Handpiece Oil

It is easy and accurate to clean and fill the dental handpiece at any moment. It has high safety. Easy to operate and easy to master. The utility model can simultaneously clean and grease 3pcs handpiece, and the rotating gear can effectively clean and grease the handpiece from different angles so as to complete cleaning and oiling

Operating volatile: 220V/110V Optional
Import pressure scope: 4-8kg/cm2
Working pressure: 2.2-2.5kg/cm2
Volume of groove cavity: 300ML
Net weight: 7KG
package size: 32*26.5*38.5



MYD-064
银汞调和机
Amalgamator



MYD-065
石膏震动调合机
Plaster rocker



MYD-066
超声波清洗机 (1.3L)
Digital ultrasonic cleaner



MYD-067
超声波清洗机 600ML
Ultrasonic cleaner



MYD-068
超声波清洗机
Ultrasonic cleaner



MYD-069
超声波清洗机 5L
Ultrasonic cleaner



MYD-060
UV- 灭菌器
UV-Sterilizer

功率: 15W
电压: 220V 50HZ
容量: 27L
消毒时间: 15MIN Outline dimension
外形尺寸: 30*23*47CM
包装尺寸: 36*28*59CM

Power: 15W
Voltage: 220V 50HZ
Capacity: 27L
Disinfection time: 15MIN Outline dimension
size: 30*23*47CM
Packing size: 36*28*59CM



MYD-061
洗片小暗箱
Dark room



MYD-070
204 技工打磨机
204 Technical micro motor



MYD-071



MYD-072
豪华四孔喷砂枪
Luxury air polisher



MYD-062
头戴式可调焦放大镜
Head mounted zoom lens set



MYD-063
头戴式可调焦放大镜
Head mounted zoom lens



A 金属箱包装
Metal box packaging



B 纸盒包装
Carton packaging

装箱单 PACKING LIST
放大镜 main glass 1set
头绳 headband 1pcs

螺丝刀 screwdriver 1pcs
清洁布 cleaning cloth 1pcs
护镜盖 protective cap 1set



MYD-073
90 门诊打磨机
90 Clinic micro motor



MYD-074
喷砂粉
Powder

MYD-075
喷砂粉
Powder

MYD-076
喷砂枪嘴
Polisher gun



MYD-077 打磨机轴承
Micro motor bearing

Specification

- 1.LED teeth whitening machine floor stand model.
2. Two groups of LED lamps, each 4pcs of upper and lower.
3. Optional output light model: A: upper light, B: below Lights, AB: both of the two groups.
4. Output power can choose from 40% to100%
5. Optional whitening time: 10 to 40 mins.
6. Dental led bleaching system and whitening light lamp accelerator teeth
- 7.wavelength: 460-490nm
- 8.Input voltage:AC110 to 240V,50/60Hz; Watt:40W.
- 9.The whiten hour:10-15minutes for one time, two or three times for one period of the treatment.



牙科高速涡轮手机

LED E-GENERATOR HIGH SPEED HANDPIECE



MYD-H01

产品描述：风光高速手机，3孔喷水，防回吸
手机接口气压：206-231Kpa(2.1-2.35kgf/cm²)
综合治疗机压力表气压：277-316Kpa(2.8-3.2kgf/cm²)
水压：198Kpa(2kg)
转速：300,000~320,000rpm
机头高度：13.02mm
机头直径：12.5mm
夹针方式：按压式
适用车针：1.59mm~1.6mmx21mm~23mm(直径x长度)
噪音：≤70dB

Product Description: LED E-Generator high speed handpiece, triple water spray, Single water spray handpiece
Interface air pressure of HP: 206-231Kpa(2.1-2.35kgf/cm²)
Air pressure at the gauge in the dental unit: 277-316Kpa(2.8-3.2kgf/cm²)
Water pressure: 198Kpa(2kg)
Speed: 300,000-320,000rpm
Head height: 13.02mm
Head diameter: 12.5mm
Clamped way: push button
Applicable burs: 1.59mm-1.6mmx21mm-23mm(diameterx Length)
Noise: ≤70dB

标准头高速手机

STANDARD PUSH BUTTON HIGH SPEED HANDPIECE



MYD-H02

产品描述：标准头高速手机，按压式、锁针式，单孔喷水
气压：180-250Kpa(1.8-2.5kgf/cm²)
水压：198Kpa(2kg)
转速：370,000~420,000rpm
机头高度：13.14mm
机头直径：Φ11.2mm
夹针方式：按压式，锁针式
适用车针：1.59mm~1.6mmx21mm~23mm(直径x长度)
噪音：≤70dB

Product Description: push button high speed handpiece, wrench type handpiece, anti-retraction system, Single water spray handpiece
Air pressure: 180-250Kpa(1.8-2.5kgf/cm²)
Water pressure: 198Kpa(2kg)
Speed: 370,000-420,000rpm
Head height: 13.14mm
Head diameter: 11.2mm
Clamped way: push button, key lock
Applicable burs: 1.59mm-1.6mmx21mm-23mm(diameterx Length)
Noise: ≤70dB



MYD-H03

大风轮高速手机（按压式、三点喷雾）
Torque big head high speed handpiece (push button)



MYD-H04

迷你头高速手机（取针式）
Mini high speed handpiece (wrench type)



MYD-H05

大风轮高速手机（取针式）
Torque high speed handpiece(wrench type)



MYD-H06

水珠款高速手机
High speed handpiece(wrench type)



MYD-H11



MYD-H12



MYD-H13



MYD-H14



MYD-H07

TORQUE TRIPLE SPARY HIGH SPEED HANDPIECE

Feature: Triple spary, anti-retraction.
Chuck type: wrench and push button
Handpiece head: torque
Air pressure: 180-250Kpa(1.8-2.5kgf/cm²)
Water pressure: 198Kpa(2kg)
Speed: 370,000-420,000 rpm
Head height: 13.14mm
Head diameter: 11.2mm
Noise: ≤70dB



MYD-H08

QUICK COUPLING HIGH SPEED HANDPIECE

Feature: Triple spary, anti-retraction, quick coupling
Chuck type: wrench and push button
Spary: triple water spary
Handpiece head: Standard
Working air pressure: 0.20MPa-0.30MPa
Rotation speed: 350,000-4000,000 R.P.M
Water pressure: 198KPa (2kg)
Noise: ≤70dB



MYD-H09

TORQUE TRIPLE SPARY HIGH SPEED HANDPIECE

Feature: Triple spary, anti-retraction.
Chuck type: push button
Handpiece head: torque
Air pressure: 180-250Kpa(1.8-2.5kgf/cm²)
Water pressure: 198Kpa(2kg)
Speed: 370,000-420,000 rpm
Head height: 13.14mm
Head diameter: 11.2mm
Noise: ≤70dB



MYD-H10

45 DEGREE HIGH SPEED HANDPEICE

Feature: Single spary, anti-retraction, 45 degree
Chuck type: wrench and push button
Spary: single water spary
Handpiece head: Standard
Working air pressure: 0.20MPa-0.30MPa
Rotation speed: 350,000-4000,000 R.P.M
Water pressure: 198KPa (2kg)
Noise: ≤70dB



MYD-H15



MYD-H16



MYD-H17



MYD-H18



MYD-H19



MYD-H20



MYD-H21



MYD-H22

Description: Low speed handpiece set 1:1 ratio, outside cooling systems
Air pressure: 245-392Kpa (2.5-4.0Kgf/cm²)
Water pressure: 198Kpa (2kg)
Rotate speed: 22,000-27,000rpm
Bur applicable: 2.334-2.355mm
Noise: <70Db



- Specifications**
- | | |
|-----------------------|------------------------------|
| Input Voltage: | AC 12-24V |
| Power : | 28W |
| Illumination: | 8000-2300Lx |
| Color Temperature | 4500K |
| Spot Size: | 80*160mm(Distance for 700mm) |
| Amount of LED tube: | 9pcs |
| Dimension of carton: | 44.5*38.5*20.5cm |
| Weight of one carton: | 2.7KG |
| Amount of one carton: | 1piece |
1. Exclusively obtain 3-len imaging technique with patent right.
 2. Square, neat light spot
 3. Illuminance within the light spot is average, light's cutting-off show very good
 4. Eight degree of adjustable illuminance with sensor, with the buzz voice indicator
 5. With lighting bands to show the illuminance
 6. 9 LED bulbs, adopt "CREE" IC made in USA
 7. Appearance design show noble and graceful, protected by patent right



- Specifications**
- | | |
|-----------------------|--------------------------------|
| Input Voltage: | AC 12-24V |
| Power: | 20W |
| Illumination: | 8000-20000Lx |
| Color Temperature | 4500+3000K(P103A)/4500K(P103B) |
| Spot Size: | 80*160mm(Distance for 700mm) |
| Amount of LED tube: | 6pcs |
| Dimension of carton: | 44.5*38.5*20.5cm |
| Weight of one carton: | 2.7KG |
| Amount of one carton: | 1piece |
1. Exclusively obtain 3-len imaging technique with patent right.
 2. Square, neat light spot
 3. Illuminance within the light spot is average, light's cutting-off show very good
 4. Two degrees of color temperature available for choice, that is 4500K and 3000K, 3000K is the yellow light is designed specially for filling teeth with resin, which is not easy to be solidified
 5. Eight degree of adjustable illuminance with sensor, with the buzz voice indicator
 6. Streamlined design with patent right protecting



- Specifications**
- | | |
|-----------------------|------------------------------|
| Input Voltage: | AC 12-24V |
| Power: | 15W |
| Illumination: | 6000-15000Lx |
| Color Temperature | 4500K |
| Spot Size: | 80*160mm(Distance for 700mm) |
| Amount of LED tube: | 4pcs |
| Dimension of carton: | 39.5*33.5*17.5cm |
| Weight of one carton: | 2.1KG |
| Amount of one carton: | 1piece |
1. Exclusively obtain 3-len imaging technique with patent right.
 2. Square, neat light spot
 3. Illuminance within the light spot is average, light's cutting-off show very good
 4. Eight degree of adjustable illuminance with sensor, with the buzz voice indicator
 5. Appearance design look compact and graceful





MYD-L07
豪华口腔灯 (带感应)
luxury operating light (with sensor)



MYD-L08
圆形口腔灯 (带感应)
round operating light (with sensor)



MYD-L09
方形口腔灯 (塑料)
square operating light (plastic)



MYD-L10
LED 口腔灯四珠
LED light (4 bulb)



MYD-L11
LED 口腔灯六珠
LED light (6 bulb)



MYD-L12
LED 口腔灯四珠
LED light (4 bulb)



MYD-L13
LED 口腔灯六珠
LED light (6 bulb)



MYD-L14
LED 口腔灯四珠
LED light (4 bulb)



MYD-L15
LED 口腔灯四珠
LED light (4 bulb)



MYD-L16
LED 口腔灯六珠
LED light (6 bulb)



MYD-L17
LED 口腔灯反射式
LED light (Reflection type)



MYD-L18
LED 口腔灯反射式
LED light (Reflection type)



MYD-L19
灯臂
Arm



MYD-L20
中柱
Post



MYD-L21
圆形反光镜
Circular mirror



MYD-L22
方形反光镜
Square mirror



MYD-L23
圆形灯罩
Round lens



MYD-L24
方形灯罩
Square lens



MYD-L25
灯泡护罩
Bulb shield



MYD-L26
口腔灯开关
Light switch



MYD-L27
口腔灯泡
Bulb



MYD-L28
带线口腔灯泡
Bulb with wire



MYD-L29
灯座
Lamp holder



MYD-L30
灯座套管
Light base holder



MYD-D01
医师座椅
Dentist stool



MYD-D02
医师座椅
Dentist stool



MYD-D03
医师座椅
Dentist stool



MYD-D04
医师座椅
Dentist stool



MYD-D05
医师座椅
Dentist stool



MYD-D06
医师座椅
Dentist stool



MYD-D07
医师座椅
Dentist stool



MYD-D08
医师座椅
Dentist stool



MYD-D09
进口脚轮
Import castor



MYD-D10
医师椅座垫
Cover cushion



MYD-D11
医师椅靠背
Backrest



MYD-D12
医师椅大圆盘
Large disk



MYD-D13
医师椅底盘
Turning chassis



MYD-D14
医师椅转盘
Turing chassis



MYD-D15
铸铝五爪盘
Metal leg



MYD-D16
尼龙五爪盘
Nylon leg



MYD-D17
黑色气弹簧
Gas spring



MYD-D18
电镀气弹簧
Gas spring



MYD-D19
静音脚轮
Foot wheel



MYD-D20
普通脚轮
Foot wheel



MYD-F01 水调节阀 (5mm)
Regulate valve (l)
MYD-F02 水调节阀 (5mm)
Regulate valve (l)
MYD-F03 水调节阀 (3mm)
Regulate valve (s)
MYD-F04 水调节阀 (3mm)
Regulate valve (s)
MYD-F05 水转换开关 (拨动)
Water exchange valve
MYD-F06 水转换开关 (拨动)
Water exchange valve



MYD-F07 总气开关
Main air switch
(3mm or 5mm)
MYD-F08 水转换开关 (推拉)
Water exchange valve
MYD-F09 水转换开关 (旋钮)
Water exchange valve
MYD-F10 两孔圆形脚踏阀
Foot pedal valve
MYD-F11 四孔圆形脚踏阀
Foot pedal valve
MYD-F12 工作状态阀
Normal valve



MYD-F13 长形四通 (3mm)
Square connector
MYD-F14 长形四通 (5mm)
Square connector
MYD-F15 吹屑阀
Blowing valve
MYD-F16 吹屑阀
Blowing valve
MYD-F17 常开挂架阀 (大)
Normal open hanging valve
MYD-F18 常闭挂架阀 (大)
Normal close hanging valve



MYD-F19 常开挂架阀
Normal open hanging valve
MYD-F20 常闭挂架阀
Normal close hanging valve
MYD-F21 单联吸吐阀
Suction valve
MYD-F22 强吸阀
Strong suction valve
MYD-F23 弱吸阀
Weak suction valve
MYD-F24 四孔方形脚踏阀
Foot pedal valve



MYD-F25 四孔方形脚踏阀
Foot pedal valve
MYD-F26 脚踏阀
Foot pedal valve
MYD-F27 两孔方形脚踏阀
Foot pedal valve
MYD-F28 洁牙机快接头 (防回吸)
Scaler quick connector
MYD-F29 洁牙机快接头
Scaler quick connector
MYD-F30 箱体多通
Integration connector



MYD-F31 箱体多通
Integration connector
MYD-F32 长形六通
Square connector
MYD-F33 副控挂架阀
Hanging valve
MYD-F34 上挂挂架阀
Up hanging valve
MYD-F35 上挂挂架阀
Up hanging valve
MYD-F36 吸吐阀
Suction Valve
MYD-F37 吸吐阀
Suction valve



MYD-F38 六角压水阀
Pressure Valve
MYD-F39 六角压水阀
Pressure Valve
MYD-F40 六角压水阀
Pressure Valve
MYD-F41 方形压水阀
Pressure Valve
MYD-F42 方形压水阀
Pressure Valve
MYD-F43 箱体集成阀
Integration valve
MYD-F44 地箱阀
Floor Box Valve



MYD-F45 四联膜片阀 (佛山)
4 In 1 valve
MYD-F46 四联膜片阀 (上海)
4 In 1 valve
MYD-F47 三联膜片阀 (佛山)
3 In 1 valve
MYD-F48 箱体集成阀
Integration valve
MYD-F49 简易减压阀
Simple pressure reducer
MYD-F50 自来水转换接头
Tap water transfer connection



MYD-F51 换向阀
Direction valve
MYD-F52 气锁开关
Air switch
MYD-F53 进气开关
Induction switch
MYD-F54 进水开关
The water switch
MYD-F55 气控电开关
Electrical valve
MYD-F56 单向阀 (5mm)
One way valve
MYD-F57 单向阀 (3mm)
One way valve




MYD-F58 四孔螺旋 (4 件套)
4 Hole connector
MYD-F59 两孔螺旋 (4 件套)
2 Hole connector
MYD-F60 四孔螺旋
4 Hole connector
MYD-F61 两孔螺旋
2 Hole connector
MYD-F62 四孔手机接头
4 Hole connector
MYD-F63 两孔手机接头
2 Hole connector
MYD-F64 气刹阀
Air brake



MYD-F65 水过滤器 (铜)
Water filter
MYD-F66 水过滤器 (塑料)
Water filter
MYD-F67 带杯空气减压阀
Air filter reducer
MYD-F68 带杯空气减压阀
Air filter reducer
MYD-F69 空气减压阀
Air reducer
MYD-F70 圆形压力表
Round meter



MYD-F71 220V 加热器
220 V heater
MYD-F72 24V80W 加热器
24V 80w heater
MYD-F73 铝合金强吸头
Strong suction
MYD-F74 铝合金弱吸头
Weak suction
MYD-F75 带过滤杯
Belt filter cup
MYD-F76 方形压力表
Square pressure meter





MYD-F144
五联挂架
Handpiece holder with frame



MYD-F145
五联挂架
Handpiece holder with frame



MYD-F146
五联大挂架
Handpiece holder



MYD-F147
五联挂架
Handpiece holder



MYD-F148
单联挂架
Single holder



MYD-F149
新款可调强吸头
Strong suction



MYD-F150
新款可调弱吸头
Weak suction



MYD-F151
可调强吸头
Strong suction



MYD-F152
可调弱吸头
Weak suction



MYD-F153
强吸头
Strong suction



MYD-F154
弱吸头
Weak suction



MYD-F155
工具盘集污杯
Dirt collector



MYD-F156
弱吸接头
Weak suction connector



MYD-F157
强弱吸转接头
Strong to weak connector



MYD-F158
工具盘拉手
Tray holder



MYD-F159
冲痰弯管
Water flush tube



MYD-F160
供水弯管
Water flush tube



MYD-F161
金属冲痰管
Water flush metal tube



MYD-F162
冲痰供水座
Water supply holder



MYD-F163
强吸嘴
Strong suction tips



MYD-F164
排水三通
Adaptor



MYD-F165
排水三通
Adaptor



MYD-F166
排水多通
Multi adaptor



MYD-F167
过滤器
Filter



MYD-F168
滤网滤盖 5 件套
Spittoon filter set



MYD-F169
滤网滤盖 4 件套
Spittoon filter set



MYD-F170
痰盂滤网 (高 / 矮款)
Spittoon filter (high or low)



MYD-F171
手机管装饰件
Handpiece tube cover



MYD-F172
过滤杯
Filtering cup



MYD-F173
弱吸集污杯
Weak suction cup



MYD-F174
塑料三通
Three way adaptor



MYD-F175
塑料三通
Three way adaptor



MYD-F176
塑料三通
Three way adaptor



MYD-F177
塑料三通
Three way adaptor



MYD-F178
塑料四通
Four way adaptor



MYD-F179
塑料四通
Four way adaptor



MYD-F180
塑料直通
Adaptor small



MYD-F181
塑料直通
Adaptor big



MYD-F182
手机管装饰套
Cell phone tube decoration



MYD-F183
快捷三通
Quick tee



MYD-F184
快捷直通
Quick and direct



MYD-F185
牙椅坐垫
Seatrest cushion



MYD-F186
牙椅靠背
Backrest cushion



MYD-F187
牙椅头枕
Pillow cushion



MYD-F188
进口光纤手机管
Fiber optic
Handpiece tubing (import)



MYD-F189
4 孔手机管 (进口)
4 hole handpiece tubing (import)



MYD-F190
光纤盒
Fiber box



MYD-F191
光纤盒 (进口)
Fiber box (import)



MYD-F192
水气管
Tube

Size	Mroll
8×5	100
6×4	200
4×2.5	200
3×2	200



MYD-F193
光纤手机管
Fiber optic
Handpiece tubing



MYD-F194
4 孔手机管带头 (抗油)
4 hole with connector



MYD-F195
4 孔手机管 (抗油)
4 hole handpiece tubing



MYD-F196
2 孔手机管带头 (抗油)
2 hole with connector



MYD-F197
2 孔手机管 (抗油)
2 hole handpiece tubing



MYD-F198
4 孔复合管 (带头)
4 hole with connector



MYD-F199
4 孔复合管
4 hole handpiece tubing



MYD-F200
4 孔手机管带接头
4 hole with connector



MYD-F201
4 孔手机管
4 hole handpiece tubing



MYD-F202
2 孔手机管带接头
2 hole with connector



MYD-F203
2 孔手机管
2 hole handpiece tubing



MYD-F204
4 孔螺旋管 (带接头)
4 hole with connector



MYD-F205
4 孔螺旋管
4 hole handpiece tubing



MYD-F206
2 孔螺旋管 (带接头)
2 hole with connector



MYD-F207
2 孔螺旋管
2 hole handpiece tubing



MYD-F208
三用枪螺旋管
syringe tubing



MYD-209
三用枪联管
syringe tubing



MYD-F210
三用枪直管
syringe tubing



MYD-F211
4 孔脚踏开关联管
4 Hole foot control tubing



MYD-F212
4 孔脚踏开关管 (硅胶)
4 Hole foot control tubing



MYD-F213
2 孔脚踏开关管
2 Hole foot control tubing



MYD-F214
强吸管
Strong suction tubing



MYD-F215
弱吸管
Weak suction tubing

